


HydroLead Analyzer



Operational Description:

Data can be downloaded via RS-232C, USB, or wireless data transmission protocols defined by the user.

Specifications:

- Automated or On-Demand Analysis
- Sample results in less than 10 minutes
- Total sampling and analysis costs less than \$1.50 per sample
- Significantly lower compliance costs
- Simple, one-touch operation
- Autocalibration and QA/QC protocols available
- Detection Limit -
- Weight
- Dimensions
- Power Requirements 110 VAC/240 VAC/12 VDC

For more information, or to inquire about other metals methods, please contact:



Environmental & Life Support Technology, Inc.

495 Uinta Way, Suite 210

Denver, CO 80230

www.elstechnology.com

cliffjolly@elstechnology.com

+1.303.840.4221